

Sreening of surface

- Switch on the two heaters.
- Hold the glasses with its` temples turned up in front of the black contrast-plate.
- Twist the glasses up to 360°.

Inspection of tension

- Place the polarisation-plate on the bearing surface.
- Hold the glasses between polarisation-plate and face-plate.

Discovering engravings

Transmitted light

- Illumination from behind
- Use contrast-plate 1.
- Use the magnifier.


Illumination from below (without picture)


- Ilumination from below
- Use contrast-plate 2 (black).


Reflexion


- Illumination from below / behind
- Hold the glasses in a way that light is directly reflected from the surface to the observer.











Inspection of centering

The scale-plate consists of the following scales:

Bridge-Scale

Helps to arrange the glasses.

Single Vision-Scale

To determine the single PD of single-vision lenses.

Bifocal-Scale

To dertermine the near PD. Nasal or lateral area can be used for measuring.

Progressive-Scale

Both of the progressive scales are spaced at an interval of 34 mm. This enables the measurement at the nasal or lateral engraving.

Height-Scale

On this scale the height of single-vision, bifocal- an progressive lenses can be measured.

Marking Zone (optional for Q-Check002)

To mark the centration point of progressive lenses.

Progressive lenses

- 1. Put contrast-plate 1 onto the transparent bearing surface.
- 2. Align glasses.
- 3. Use the magnifier.
- The position of the engraving on one of the progressive-scale = PD/2
- The height read off + interval value stipulated by the lens manufacturer engraving ↔ centering cross = height (regarding frame edge)

Bifocal-lenses

The near PD

- 1. Put the polarisation-plate onto the bearing surface.
- 2. Adjust the glasses in such a way that the segment is in centre of the bifocal-scale = value 1.
- Move glasses until they are adjusted with the bridge in centre. Read off the value of the same segment = value 2.
- 4. Difference = value 1 value 2
- 5. Near PD = 30 difference

Height

Read the height-scale.

Example

Value 1 = 34,5Value 2 = 33Difference = 34,5 - 33 = 1,5Near PD = 30 - 1,5 = 28,5Height = 19