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Multifunctional instrument for quality inspection of lenses

Standards of quality and precision of lenses have risen over the last decades. Particularly video lens centering systems, CNC-Technology and individual progressive lenses have contributed to the fact that the customer nowadays gets high quality and precise lenses from opticians. Innovation in combination with an improvement of quality only makes sense, if it can be inspected. Unfortunately the quality inspection of the improved quality lenses has been neglected in opticians' workshops. Over the last decades inspections were made with a PD ruler – a method which is in these days not suitable for „high tech-lenses“.

In particular progressive glasses that come along with high production requirements could not be inspected regarding their centering. Instead of a direct inspection the centering could only be tested by using the engravings and a suitable pattern to mark the centration point and then finally measure the centering – a time consuming procedure. The multi-level process described is pedestrian and because of its many steps makes a lot of mistakes possible. Every step of the process is bound to inaccuracies. An Accumulation of these small mistakes leads to a rather inexact result as all of the single inaccuracies interfere with each other.

The optical effects of the lens that occur while looking through it and the parallax errors will reinforce the „wrong“ result. All of this makes the multi-level inspection of centering rather unstable.

At opticians' suggestions we have dealt with the problem and can present Q-CHECK an instrument for the multifunctional quality inspection of lenses. Our main aim was to develop a practical device that enables the user to inspect quickly and precisely. Being a compact all-rounder Q-CHECK can ideally be used in the optician's workshop for a final inspection. Furthermore Q-CHECK eases an inspection of centering and coating defects in the process of dealing with complaints. In the service department the optician can convince his customers of how precisely the lenses are manufactured.

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The multifunctional Q-CHECK ensures an extensive quality check of lenses:

Inspection of centering:

Q-CHECK allows the measuring of single-vision and multifocal lenses. A direct inspection is possible with progressive lenses – a time consuming and incorrect marking is unnecessary. The inspection is performed with the engraving that is stable and unchanging at a standardized position. In addition to that the centration point can be marked clearly and exactly for an inspection later on together with the customer (Q-CHECKplus) – without marking the engraving and without any of the mistakes described above. Looking for a suitable and compatible pattern is unnecessary as well as the Q-CHECK scale is compatible with progressive addition lenses of several manufacturers! Also a progressive addition lens that is a little twisted within the customer's spectacle lens can be detected.

Discovering engravings:

Because of a special contrast scheme and the magnifier's increasing impact, discovering engravings is very easy. In addition to those features Q-CHECK offers several methods of discovering engravings.

Inspection of tension:

Due to big, high-quality polarization foils an inspection of both spectacle lenses is possible at the same time – quick and effective.

Surface screening:

The surface screening enables the optician to inspect lenses and detect material and surface defects on equal terms and under repeatable conditions – needless to say essential for providing high quality.

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Safety Instructions

Intended Use

This device has been developed for the quality inspection of lenses regarding their centering, surface defects and of tension according to RAL-RG 915. It is not to be used for any other reason.

General safety

- Carefully read the safety information and the operating instructions before connecting up the unit. This will ensure that you use all functions safely and reliably.
- Retain these operating instructions.
- Ensure that these instructions are passed on should the unit be transferred.
- Packaging materials do not constitute a toy. Never allow children to play with plastic bags as there is a risk of suffocation.
- Never attempt to open or repair the unit yourself. Only permit authorized personnel to carry out technical repairs.

Ambient and installation conditions

- Make sure that the unit is installed on a solid surface and cannot tip over.
- Never operate the unit in an environment in which water or liquid could penetrate into the unit. There is a risk of electrocution should water or liquid penetrate into the unit.
- Make sure that the cable is laid in a way that ensures that no one trips over it.
- The slots and openings on the unit are provided for ventilation. Never close these openings as the unit could otherwise overheat.

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Precautions against electric and magnetic dangers

- Make sure that the device is connected with a bus bar that provides the same voltage and frequency as stated on the type plate and that it is grounded according to the local requirements.
- The device is only to be run with a voltage of 230 V/50 Hz.
- The device has a charging rate of 0.11 A.
- The device has (on full power) a power input of 25 W.
- To pull the plug, always pull the plug, not the cable.
- For a better protection against the dangers of electric current a special safety switch should be used if an installation for that purpose is missing.
- The device contains magnets (as in the magnifier). Please be aware that magnets may delete data media and influence on electronical and mechanical components such as e.g. pacemakers. Do pay attention to the given security distances provided by the manufacturers of these products. In the case of direct contact with magnetic material and magnetic systems allergic reactions may arise (e.g. with ceramic and metallic materials as well as zinc, nickel and plastics).
- Keep the magnets away from children. Accidentally swallowed magnets can result in serious medical problems.

Cleaning and maintenance

- Before any cleaning or maintenance, switch the device off and disconnect it from the bus bar by pulling the plug out of the socket.
- Never switch on the device during cleaning or maintenance.

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1 x Q-CHECKplus Art. Nr. 001.393

In case of interferences

- Pull the plug immediately if cable or device are damaged.
- Do not induct any objects into the device. Conductor rails might be touched.
- If liquids or foreign substances get into the device, pull the plug immediately. The device should then be checked before it is put into operation again.
- If a proper function is not possible even if the function is according to all articles of the operation manual, the device has to be disconnected.

Besides these advices in the operation manual all local rules for accident prevention have to be considered and followed.

Scope of supply and description of appliance

Q•CHECK allows a precise, effective and quick quality inspection of mounted glasses according to RAL-RG 915. The use of various contrasts and illuminating procedures makes 4 different functions possible. Q•CHECK offers the inspection of centering of single-vision-, bifocal- and progressive lenses, eases the locating of glass-engravings, provides a tension inspection and enables the user to make a surface screening.

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1 x Lupe Q-CHECKplus Art. Nr. 000.775

Q•CHECKplus

- 1 Transparent back
- 2 Face-plate with scale
- 3 Transparent bearing surface
- 4 On-off switch
- 5 Contrast-plate (black, fixed)
- 6 Heater (surface)
- 7 Two parallel arranged heaters (2700 K)
- 8 Magnifier

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1 x Contrast-plate 1 Q-CHECK Art. No. 000.661
(for the purpose of discovering engravings - transmitted light)

1 x Plaque de contraste 1 Q-CHECK Art. No. 000.661
(pour la reconnaissance de gravure - en transparence)

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1 x Contrast-plate 2 Q-CHECK Art. No. 000.662
(for the purpose of discovering engravings - illumination from below)

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1 x Clip Q-CHECKplus Art. No. 0001.235

Please make sure your set is complete. If any parts are missing please contact your sales partner or LUMOS (info@lumos-optics.de, Tel.: +49 (2103) 7948-626).

Set-up and required space

Please always follow these instructions:

- The device must be placed on a solid and stable sub-surface.
- Do not run the device in surroundings where water or liquids might enter. In case of ingress of water there is certain risk of electric shock.
- Make sure that the cable is arranged in a way that no one might stumble.

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Functional scope and operation

By using Q•CHECK the following inspections can be performed:

- Discovering engravings
- Inspection of centering of
 - glasses with single-vision lenses (mark with lensmeter beforehand)
 - glasses featuring progressive lenses
Marking the centration point of progressive lenses (optional with Q•CHECKplus)
 - glasses with bifocal- und trifocal-lenses
 - edged lenses
- Inspection of tension
- Surface screening regarding material and surface defects
- Checking the arrangement of polarized lenses in mounted glasses

Please be aware that the illuminants used in this device take a certain amount of time (approx. 5 min.) before operating at full capacity.

Discovering engravings

Depending on the type of glass, production process, surface coating and manufacturer engravings can be discovered and made visible under various conditions.

Q•CHECK therefore provides three different ways of discovering engravings.

- Transmitted Light in front of bright-dark-contrast
- Illumination from below in front of a black Contrast-plate
- Discovery in reflexion

The magnifier that is held magnetically at the black lower casing facilitates the discovery of the engraving and reading off its precise location. In addition to the magnifier provided a stand magnifier with a bigger range of vision may be used. A trial lens can be inserted into the lens holder so that a presbyop can recognize the scale. Not only older but also younger users will find the addition comfortable to read off the scale.

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Transmitted Light: Illumination from behind through the lenses; use Contrast-plate 1

This technique works with relatively „flat“ engraving structures which can be found in several synthetic lenses as used by major manufacturers. It allows direct reading off the location of the engraving and according to that the value of centering.

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Illumination from below: Use Contrast-plate 2

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In case of „steep“ engravings this technique should be chosen. It allows the direct reading off and pointing out the centering values by using the scale as well.

Reflexion: Hold the glasses in a way that light is directly reflected from the surface to the observer.

Engravings can be successfully discovered and found as the reflection differs from the surrounding to the engraving. As a matter of principle this technique demands the marking out of the engraving with a pen.

The magnifier is fixed magnetically to the lower casing to ease the discovery of engravings.

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Inspection of centering

The inspection of centering of all types of lenses is accomplished by using the scale. The scale is arranged in a way that the values of centering can even on multi-focal lenses be read off immediately and regarding that is the most important element of the inspection of centering.

The scale-plate consists of the following scales:

Bridge-scale: Helps to arrange the glasses in Q•CHECK.

Single vision-scale: To determine the Single PD of Single-vision lenses.

Bifocal-scale: To determine the near PD. The centers of the Bifocal-Scales are shifted outwards 30mm, as this is a common near PD value on which the measuring of the near PD (explained later on) depends. Nasal or lateral area can be used for measuring. Both symbols in shape of a segment define the region where the segment is to be aligned.

Progressive-scale: Both of the progressive scales are spaced at an interval of 34 mm. This marks the distance between both of the engravings within the progressive lens. (manufactured according to EN ISO 8980-2). This enables the measurement at the nasal or lateral engraving. The upper circles „nasal“ and „lateral“ symbolize this scale.

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
Height-scale: On the Height-scale the height of single-vision, bifocal- and progressive lenses can be measured.

Marking zone: To mark the centration point of progressive lenses.

Even though the scale appears a bit unfamiliar at first one gets used to its sections. Dotted lines facilitate the discovering of 32 as „middle PD“, the heights in intervals of 10mm, as well as alignment of the segment. The broader vertical lines are arranged in intervals of 2mm, the finer lines are arranged in intervals of 1mm. The height-lines are spaced in intervals of 2mm.

Q-CHECK enables the user to inspect the centering of manufactured glasses with single-vision-, bifocal- and progressive lenses.

Positioning and arrangement of glasses

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
Adjustment of the glasses

The glasses should be adjusted with temples spread on the bearing surface behind the face-plate (press against it). The lower frames must bear on the surface and glasses and face-plate must be in direct contact to avoid mistakes during the measuring process. To attach the glasses to the scale the clip provided can be used. It should be affixed magnetically to the lower illumination plate. By twisting the clip the glasses can be moved horizontally. When working with thin bridges the clip can be moved closer to the scale.

To perform a quick inspection the glasses can also be held manually.

Align glasses in center with upper and lower scale.

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
Contrast plates: All contrast plates provided facilitate a discovering of engravings in progressive lenses. These are preferably used during the following inspections:

Contrast plate 1: Discovering engravings (transmitted light)

Contrast plate 2: Discovering engravings (illumination from below)

Polarization plate: This plate makes an inspection of tension and the inspection of bifocal lenses possible.

Changing of the plates

All plates can be used depending on the contrast procedure and stick magnetically to the illuminated rear area. Thus the plates can easily be removed and changed. The plate-holder at the aft casing can be used to store the plates which are not being used at that point of time.

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Reading off the scale

To work parallax-free the scale should be read off in compliance with the measuring method of Viktorin!

The reflection of the Q•CHECK sticker must be vertical in the center of the magnifier. It must also be in coincidence with the read off PD-value.

With Q•CHECKplus the magnifier can be moved horizontally by moving the slider.

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

Inspection of centering

Lenses

All lenses must possess marked points of reference (focimeter).

1. No special plate is necessary.
2. Align glasses
3. At the „single vision scale“ left and right PD values and at the Height-scale the centering height (according to the lower frame) can be read off. To calculate the height according to the top of the bevel, the remaining material thickness from the lower frame in the groove must be calculated and subtracted beforehand.

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Progressive lenses

As a matter of principle two different forms of inspection must be pointed out:

1. Inspection of centering via its „pressure“

Should the pressure of the inspected progressive lenses be adequate, then inspection can proceed according to inspection requirements for single vision lenses.

2. Inspection of centering via engravings

In many cases such an inspection will prove as impossible or not reasonable as the pressure is no longer present. Also an inspection via engravings will lead to a result of higher precision.

It is a precondition for the inspection via engravings that all engravings are to a certain extent visible which might not be the case with some of the older lenses. In that case the engravings must be marked beforehand (see picture).

More detailed explanations on how engravings can be detected with Q•CHECK can be found in chapter „Discovering engravings“.

A check can be made by using the progressive scale which consists of a nasal and lateral section. The progressive scales are spaced at an interval of 34 mm. This value is compliant with the interval of both engravings in a progressive lens (manufactured according to EN ISO 8980-2). Thus it is possible to measure at the nasal and lateral engraving.

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Glasses with progressive lenses can if they don't feature any markings be inspected via engravings.

1. Put contrast-plate 1 onto the transparent bearing surface. (Alternatively: contrast-plate 2)
2. For a better discovering of the engravings use the magnifier. It can be fixed magnetically to the black, lower case. The way he holds his hand (see picture) makes it easy to measure as the other hand can be used to adjust the magnifier.
3. Align glasses. Push glasses at the bridge against the scale. If necessary use the clip.

4. The nasal and lateral engraving can be read off. Make sure that both of them are on the same height if the lens has just been edged before. The center of the engraving on the nasal and lateral scale complies with the single-PD.

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




5. To read off the centering height add the vertical interval between engraving and centering-cross which is depending on the value stipulated by the lens manufacturer of.

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


6. All heights refer to the lower frame! To calculate the height according to the top of the bevel, the remaining material thickness from the lower frame in the groove must be calculated and subtracted beforehand.

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

Marking the centration-point of progressive lenses (optional with Q•CHECKplus)

Via the marking section at the upper frame of the scale the centration point can easily and securely be marked. This can be very useful to verify the centering once more before handing out the glasses to the customer or to determine the centering in the edging process. Align the lens and both of the engravings with the scale. The engravings must be on the horizontal 32 line of the progressive scale. The height depends on the engraving centering-cross by the lens manufacturer (see chapter above). Nasal and lateral engravings must aligned with that value. Throughout the hole the centration-point can be marked with thin pen.

Example: Vertical interval engraving - centration-cross depending on manufacturer = 4mm

1. Adapt nasal and lateral engravings to the height
2. Mark out with a thin pen throughout the hole

(We suggest you use Staedtler lumocolor permanent special thin, which has a good adhesion on even glazed „clean-coat“ surfaces, or Edding 780 white, which creates an ideal contrast to the dark background with its white color. It can also be easily removed with benzine!)

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The way the hand is held in the picture on the left eases the marking:

Hold the glasses at the upper and lower frame. That way it can easily be moved and adjusted. Push glasses slightly against the scale with one hand and use the other hand to remove the magnifier and mark by inserting a pen through the hole.

Marking on the right side: left hand holds glasses at the right lens as right hand is marking

Marking on the left side: right hand holds glasses at the left lens as left hand is marking

Additionally the glasses should be fixed at the bridge by using the clip. By twisting the clip the glasses can be moved horizontally.

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
Bifocal-lenses

Measuring the near PD:

1. Put the polarization-plate onto the bearing surface.
2. Adjust the glasses in such a way that the segment is in center of the bifocal-scale. Read off the value of the segment at the nasal and lateral bifocal-scale = value 1 (when working with rounded segments the outermost frame of the segments must be used. If necessary remove polarization-plate!)
3. Move glasses until they are adjusted with the bridge in center. Read off the value of the same segment = value 2
4. Difference = value 1 – value 2
5. Near PD = 30 – difference

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
Difference = 34.5 – 33 = 1.5

Near PD = 30 - 1.5 = 28.5

As the center of bifocal-glasses-segments, to which the near PD refers, does not show any markings as intended by the manufacturer, a precise measuring using two values must be made via shift.

Trifocal-lenses can be measured the same way.

Measuring the segment-height:

The segment-height (referring to the lower frame) can be read off according to the height-lines at the uppermost angle of the segment. As you can see in the example above the height is 19. (To calculate the height according to the top of the bevel, the remaining material thickness from the lower frame in the groove must be calculated and subtracted beforehand.)

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Edged lenses

The height of edged lenses can also be measured before the actual mounting takes place. The lens then has to be placed directly and adjusted and the height that is read off refers to the Boxing system.

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Inspection of tension

Q•CHECK allows the inspection of tension of both lenses at the same time. To inspect mounted glasses according to their tension the polarization-plate must be placed on the bearing surface. Hold the glasses between polarization-plate and face-plate.

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
Then the tension of the edged lenses can immediately be evaluated.

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Surface screening and error detection

During the surface screening material defects, coating defects, dirt and scratches can be detected. Therefore the two heaters should be switched on. Hold the glasses with its temples turned up in front of the black contrast-plate. By twisting the glasses up to 360° and viewing under a flat angle surface defects and errors can be seen. It is necessary to twist the glasses in order to detect some errors that can only be seen if the light hits them in certain angles of incidence. Uncut lenses can of course be inspected the same way.

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Inspecting the arrangement of polarizing lenses

Hold glasses with polarized lenses in front of the face-plate or behind it (no plate on bearing surface). If the glasses have been edged correctly it is possible to look through them and impossible to look through them if the glasses are turned vertically.

Further information on operation

Tips for the use of Q-CHECK:

- Using the plates to achieve contrast enhancement is dependent on the user's subjective sense. A plate that is not recommended might also possibly yield an acceptable result. This depends on the type of engraving and the lenses' thickness.
- The illumination can be switched optionally to provide ideal illuminating conditions for all glasses. Switching on both heaters does not normally disturb the user.
- If problems arise when the centering is inspected via engravings the following steps can ease the discovering of engravings:
 - Draw a circle on the lens surrounding the engraving. This makes it easier for the user to locate it and improves orientation.
 - Try out several contrast-plates.
 - In case of using contrast-plate 1 finding the ideal contrast-area can be eased by moving the head slightly vertically.
 - Mark the engraving in its center and then read off (only 2nd choice solution!)
 - Insert an addition (e.g. +1.0) to the magnifier's lens holder. That way the accommodation onto the background can be avoided and the scale can easily be recognized.

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
Replacing bulbs

The bulbs used in the Q•CHECK are extremely durable. However, should it be necessary to replace a bulb, it should only be done by a qualified electrician! The following instructions are therefore intended for use by electricians.

Disconnect the unit from the mains power supply by disconnecting the plug before commencing the work described below!

Place the Q•CHECK on its left side to replace the bulbs.

Lower bulbs

Loosen the 4 screws on the base plate to replace the lower bulbs.

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Carefully pull out the underside of the unit to the rear.

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Turn the bulb only around a quarter of a revolution counter-clockwise and remove.

Only insert new bulbs and reassemble the unit, taking care not to squeeze the cables.

Care! Only use energy-saving bulbs type MM14132 (9 watt)!

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
Rear bulbs

Place the unit on its side and loosen the four screws on the rear plate.

Carefully lift up the back of the unit.

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Turn the bulb only a little counter-clockwise and remove it. Only insert new bulbs and reassemble the unit.

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Care! Only use energy-saving bulbs type MM 71002!

Additional maintenance and cleaning work

The Q-CHECK works largely maintenance-free.

Cleaning

Surfaces can be cleaned with a soft cloth and plastic cleaner. Never use abrasive or aggressive products, such as glass cleaner or acetone.

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LUMOS GmbH